IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR Embrace USER'S GUIDE BEFORE USING THE Embrace TEST STRIPS TO TEST YOUR BLOOD SUGAR.

Please call 1-877-979-5454 for assistance if you have any questions. We are available 24 hours a day, 7 days a week; or visit our website at www.omnishealth.com.

Intended Use

The Embrace No Code Blood Glucose Test Strips are to be used with the Embrace No Code Blood Glucose Meter to quantitatively measure glucose in capillary whole blood. The Embrace System is plasma-calibrated for easy comparison to lab result. The Embrace Blood Glucose Monitoring System is intended for self-testing by persons with diabetes and by health care professionals.

Introduction

The biosensor-based Embrace Blood Glucose Test Strips employs the enzyme Glucose Oxidase that is specific for blood glucose measurement from 20 to 600 mg/dL. The arrow printed on the test strip indicates the direction of insertion. When blood from a puncture site touches the tip of the test strip; capillary action leads the blood into the reaction zone. The enzyme reacts with glucose in blood and produces an electrical signal. This signal is measured by the Embrace Meter. Test result is verbally announced and/or displayed in 6 seconds.

Reagent Composition

Each cm² of test strip contains the following ingredients in the approximate concentrations listed below:

- Glucose Oxidase (A. Niger) 7.6%
- Electron Shuttle 53.3%
- Non- Reactive Ingredients 39.1%

Warnings and Precautions

The Embrace Blood Glucose Test Strips are for use outside the body (IN VITRO diagnostic use).

- Do not use test strips after their expiration date. The expiration date is printed on the carton and on the label of the test strip bottle.
- Do not use strips that are wet or damaged.
- Do not re-use the strips.
- If your test result is below 50 mg/dL or above 250 mg/dL; perform the Control Solution Test to make sure your system is working properly and then repeat testing once more using blood from your fingertip. If the result is still the same; contact your healthcare professional immediately.
- If you have symptoms that are inconsistent with your test results and you have eliminated common errors as described in the user's guide; contact your healthcare professional immediately.
- Never make significant changes in your diabetes treatment program or ignore symptoms without consulting your physician.
- Use the Embrace Blood Glucose Monitoring System within temperature range of 50°F to 104°F.

Storage and Handling

1. Always close the vial cap tightly after removing a test strip to avoid moisture and direct sunlight.
2. Store the test strips at temperatures between 39.2°F to 86°F.
3. Store out of direct sunlight; do not refrigerate or freeze.
4. Unopened test strips are stable until the expiration date printed on the bottle when stored properly.

5. The month printed refers to the END of that month.
6. Use within 3 months after first opening.
7. Do not handle the test strips with wet or dirty hands.

Sample Collection and Preparation

Fresh finger stick capillary whole blood maybe collected into heparin or sodium EDTA test tubes and should be used within 15 minutes. Do not test on plasma or serum samples. Do not use fluoride-based preservatives.

Test Procedure

Refer to "Testing Your Blood Glucose" section in the User's Guide of the Embrace No Code Blood Glucose Meter.

Alternate Site Testing

The Embrace No Code Blood Glucose Monitoring System can also test for blood glucose from palm and forearm (alternate site testing, or AST). Physiological differences between the finger or palm and forearm may result in different blood glucose readings (1). You should not test on palm and forearm and use finger tip testing instead when:

- You have hypoglycemic unawareness
- You suspect you blood glucose level is changing rapidly because of exercise or after rapid-acting insulin injection or it has been less than two hours after a meal.
- You will be driving a car or participating in risky activities such as operating machinery. Consult your healthcare professional to decide if alternate site testing is right for you.

Expected Values

Blood glucose value depends on many factors. The ideal range for adult with diabetes is:

- 74-106 mg/dL before meals (2), and
- Less than 140 mg/dL two hours after meal (3). It is important to consult with your physician or healthcare professional to determine an appropriate target range for you.

Quality Control

Run a control solution test by following the instructions detailed in User's Guide when:

- Your test results do not agree with how you feel.
- At least once per week to verify that the meter and test strips are working properly together.
- Your test strips were stored at temperatures and humidity outside the specified storage conditions.
- Every time you open a new bottle of test strips.
- Practice your testing technique.
- You drop your meter.

Use only Embrace Control Solution. Your test results should fall within the control range printed on the test strip bottle. If test results falls outside this range; repeat the test. Results that fall outside the range maybe caused by:

- Error in performing the control solution test
- Expired or contaminated control solution
- Test Strip is damaged
- Meter malfunction If the result continues to fall out side of the printed range; please call our authorized representative by using the information printed at the end of this insert or on the box for assistance.

Problem Solving

1. Confirm the Test Strips are within the expiration date.
2. Make sure the blood completely fills the reaction zone. "Err 5" will be verbally announced or display on the Embrace Meter if there is insufficient blood. DO NOT add a second drop of blood; discard the test strip from the meter and retest using a new Test Strip.
3. Check the system by performing the Control Solutions Test.
4. Refer to the "Solving Problems" section in the User Guide for additional information. Please call our authorized representative by using the information printed at the end of this insert or on the box for assistance or your healthcare professional.

Limitations

1. DO NOT use plasma or serum samples.
2. DO NOT test on neonatal samples (new born).
3. DO NOT test on venous or arterial samples.
4. Altitude up to 2,000 meters will not affect test results.
5. Hematocrit range: 30% to 55%.
6. Triglycerides: 260 to 360 mg/dL has no significant effect on test results (normal range 36 to 165 mg/dL).
7. Icodextrin and its metabolites (maltose, maltotriose and maltotetraose) do not significantly affect test results.
8. The following will not affect test results in expected blood concentrations: ascorbic acid (Vitamin C), uric acid, or methyl-dopa.
9. Therapeutic concentration of dopamine or Ldopa may affect results. Individuals taking these drugs should not use this system to test for blood sugar.
10. Higher than therapeutic concentrations of acetaminophen, or glibenclamide may affect test results. Individuals taking these drugs should not use this system to test for blood sugar.

Performance Evaluation Data

Accuracy: Total of 158 fingertip blood samples were collected for system accuracy evaluation, 17 samples were from glucose level below 75mg/dL and samples were from glucose level above or equal to 75mg/dL. Total of 152 fresh blood samples were collected from palm and forearm for system accuracy evaluation of alternate site testing, 11 samples were from glucose level below 75mg/dL and 141 samples were from glucose level above or equal to 75mg/dL. The system accuracy evaluation of fingertip and alternate site samples were tested by healthcare professional, and compared to YSI 2300 Blood Glucose Analyzer.

Table 1: Accuracy for blood glucose level <75mg/dL

Sample Source	Within ± 5mg/dL	Within ± 10mg/dL	Within ± 15mg/dL
Capillary	7/17 (41%)	15/17 (88%)	17/17 (100%)
Palm	8/11 (73%)	11/11 (100%)	11/11 (100%)
Forearm	7/11 (64%)	10/11 (91%)	11/11 (100%)

Table 2: Accuracy for blood glucose level ≥75mg/dL

Sample Source	Within ± 5%	Within ± 10%	Within ± 15%	Within ± 20%
Capillary	78/141 (55%)	121/141 (86%)	139/141 (99%)	141/141 (100%)
Palm	81/141 (57%)	116/141 (82%)	136/141 (96%)	139/141 (99%)
Forearm	71/141 (50%)	123/141 (87%)	134/141 (95%)	141/141 (100%)

The result showed the Embrace System compares well with the YSI 2300 Glucose Analyzer and meets the minimum acceptable performance criteria defined in ISO 15197.

Precision: Three lots of Embrace No Code Test Strips were used for a within-run precision study. Venous blood in herapin tubes was spiked to 5 concentrations. Blood glucose readings were recorded for 10 days resulting in 100 data points for each concentration; as shown in the following tables:

Reagent Lot FS 001A

Number of Readings:	100	100	100	100	100
Average Reading(mg/dL):	46	113	192	303	429
S.D. (mg/dL):	2.1	3.9	6.9	9.7	13
C.V%:	---	3.4	3.6	3.2	3

Reagent Lot FS 001B

Number of Readings:	100	100	100	100	100
Average Reading(mg/dL):	46	111	189	301	422
S.D. (mg/dL):	2.7	4.3	6.3	10.8	15.4
C.V%:	---	3.9	3.3	3.6	3.6

Reagent Lot FS 001C

Number of Readings:	100	100	100	100	100
Average Reading(mg/dL):	46	109	185	296	418
S.D. (mg/dL):	3.0	4.6	7.4	11.9	13.5
C.V%:	---	4.2	4.0	4.0	3.2

References:

1. Jungheim K., Koschinsky T.: "Risky delay of hypoglycemia detection by glucose monitoring at the arm." *Diabetes Care* 2001;24(7):1303-04.
2. Stedman, Thomas Lathrop. *Stedman's Medical Dictionary, 27th Edition, 1999, pg. 2082*
3. American Diabetes Association, "Clinical Practice Recommendations 2003." *Diabetes Care, Vol 26, Supplement 1, pg. S22.*

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