

Embrace[®] EVO[™]

Compact Blood Glucose Meter



Small & Simple!
Pictured Actual Size



Easy to Use

- No coding required
- Fast – results in six seconds
- Compact design with large display
- Easy to read menu icons
- Easy to handle test strips

Test with Less Pain


- Small sample size
- Alternate site testing – fingers, palm or forearm

Manage Your Results


- 300-test memory with time/date stamp
- 7-, 14-, and 30-day averaging

EmbraceEVO[™] Test Strips

Instant Three-Step Verification ensures accurate results every time.




Confirms test strip positioned correctly



Pinpoint technology detects sample



Fill detection verifies sample is adequate



NEW! ALL-IN-ONE STARTER KITS

INCLUDES
50
TEST STRIPS

This special kit includes
50 prepackaged test strips!
Call **1-877-979-5454** for
more information!



An **Evolution** in **Affordable**, Quality Meter Technology.

EmbraceEVO[™] is a registered trademark of Omnis Health.
Copyright 2014.

Embrace[®] EVO[™] Features and Benefits

Test Strip Port

Conveniently located at the top of the meter, the test port is easy to feel for individuals with low vision. Insert a test strip and the meter will turn on automatically.

Meter Buttons

These buttons allow you to select different meter options, such as setting the time and date, programming alarms or recalling the memory.

Data Port

Connect an Omnis Health data cable here to upload your results to your PC. (Cable purchased separately).

Serial Number & Contact Information

You are important to us. We're available Monday–Friday, 8 AM–8 PM (EST). Call us toll-free at 1-877-979-5454 or visit www.embraceevo.com.

Battery Compartment


Easy to find and install, the Embrace EVO[™] uses one CR 2032 3V Lithium coin cell battery.



Meter not shown actual size.

Trust is earned.

The EmbraceEVO Meter exceeded expectations in clinical studies, achieving high ratings for accuracy, dependability and ease-of-use.*



*Studies are conducted in accordance with industry standards and are available for download at www.embraceevo.com



505 Airpark Center Dr.
Nashville, TN 37217
www.embraceevo.com

[f /embracemeter](https://www.facebook.com/embracemeter) [t @embracemeter](https://twitter.com/embracemeter)

An **Evolution** in **Affordable**, Quality Meter Technology.

MKSS00002-20140128.rev1
EmbraceEVO[™] is a registered trademark of Omnis Health.
Copyright 2014.