



Embrace®

Talking Blood Glucose Meter

Easy to Use

- No coding required
- Optional English/Spanish talking feature
- Fast – results in six seconds
- Sleek design with large display
- Easy to read menu icons

Test with Less Pain

- Small sample size
- Alternate site testing – fingers, palm or forearm


Manage Your Results

- 300-test memory with time/date stamp
- 7-, 14- and 30-day averaging
- PC downloadable results




Embrace® Test Strips


Instant Three-Step Verification ensures accurate results every time.




Confirms test strip positioned correctly



Pinpoint technology detects sample



Fill detection verifies sample is adequate



EASY. SAFE. ACCURATE.

Download test results to your computer for trending.
Software available at www.embracemeter.com.
Call **1-877-979-5454** for more information!



The **One** Choice for Talking, **Reliable** Meter Technology.

Embrace® is a registered trademark of Omnis Health.
Copyright 2014.

Embrace® Features and Benefits

Test Strip Port

Conveniently located at the top of the meter, the test port is easy to feel for individuals with low vision.

Meter Buttons

These buttons allow you to select different meter options, such as setting the time and date, programming alarms or recalling the memory.

Large Display

Whether you have low vision or not, the Embrace's large display allows easy viewing of your blood glucose results and meter functions.

Rubber Grips

The rubber grips on the side and bottom of the meter will help keep the meter from moving while testing.

Serial Number & Contact Information

You are important to us. We're available Monday-Friday, 8 AM-8 PM (EST). Call us toll-free at 1-877-979-5454 or visit www.embracemeter.com.

Data Port


Connect an Omnis Health data cable here to upload your results to your PC. (Cable purchased separately).

Battery Compartment

Easy to find and install, the Embrace uses standard AAA batteries.

Speaker


The optional talking feature of the Embrace gives results that you can hear as well as see.



Meter not shown actual size.

Trust is earned.

The Embrace Meter exceeded expectations in clinical studies, achieving high ratings for accuracy, dependability and ease-of-use.*



*Studies are conducted in accordance with industry standards and are available for download at www.embracemeter.com



505 Airpark Center Dr.
Nashville, TN 37217
www.embracemeter.com

[/embracemeter](https://www.facebook.com/embracemeter) [@embracemeter](https://twitter.com/embracemeter)

The **One** Choice for Talking, **Reliable** Meter Technology.

MKSS00001-20140128.rev1
Embrace® is a registered trademark of Omnis Health.
Copyright 2014.