

Embrace[®] PRO[™]

Professional Blood Glucose Meter



Easy to Use

- Test strip release button
- Auto coding
- Fast – results in five seconds
- Approved for multi-patient use

Test with Less Pain

- Small sample size – 0.5 μ l
- Alternate site testing – fingers, palm or forearm


Manage Your Results

- 500-test memory with time/date stamp
- Event Recording Mode for more effective data management




EmbracePRO[™] Test Strips

Accurate results within five seconds.



Side sampling: shorter conductive path increases accuracy



Gold electrode increases reliability and accuracy



NEW! PROFESSIONAL BLOOD GLUCOSE METER

Approved for multi-patient use, this meter features a **quick test strip release** for more hygienic strip disposal. Call **1-877-979-5454** for more information!



The **Professional's** choice for **Reliable** Meter Technology.

EmbracePRO[™] is a registered trademark of Omnis Health. Copyright 2014.

Embrace[®] PRO[™] Features and Benefits

Data Port
 Connect an Omnis Health data cable here to upload your results to your PC.
 (Cable purchased separately).

Meter Buttons
 These buttons allow you to select different meter options, such as setting the time and date, programming alarms or recalling the memory.

Test Strip Port
 Located at the bottom of the meter, the test port is ideal for hygienic test strip disposal.



Battery Compartment
 Easy to find and install, the Embrace PRO[™] uses one CR 2032 3V Lithium coin cell battery.


Serial Number & Contact Information
 You are important to us. We're available Monday-Friday, 8 AM-8 PM (EST). Call us toll-free at 1-877-979-5454 or visit www.embracepro.com.

Test Strip Release Button
 This feature allows for safe and hygienic disposal of a test strip after use.

Meter not shown actual size.

Trust is earned.

The EmbracePRO Meter exceeded expectations in clinical studies, achieving high ratings for accuracy, dependability and ease-of-use.*



*Studies are conducted in accordance with industry standards and are available for download at www.embracepro.com



505 Airpark Center Dr.
 Nashville, TN 37217
www.embracepro.com

[f /embracemeter](https://www.facebook.com/embracemeter) [t @embracemeter](https://twitter.com/embracemeter)

The **Professional's** choice for **Reliable** Meter Technology.

MKSS00003-20140128.rev1
 EmbracePRO[™] is a registered trademark of Omnis Health.
 Copyright 2014.